sempermed[®] (S) nitrile skin²

NITRILE SKIN² DISCOVER YOUR 2ND SKIN

PROTECTION FROM INFECTIONS – virus-proof for users and patiens¹⁾ **EXCELLENT SENSE OF TOUCH** – through optimised wall thickness **HIGH WEARING COMFORT** – through the innovative, soft material



YOUR NEEDS UNDER CONTROL

In addition to protection, wearing comfort and skin-friendliness, a good sense of touch is the most important criterion when choosing the right examination glove. The new Sempermed[®] Nitrile Skin² combines these properties in a unique, exceptionally light glove.

The soft material developed by our European Research Centre allows for maximum sense of touch and comfort due to optimised wall thickness, and at the same time maximum protection from infections. **Discover your 2nd skin!**

- Protects users and patients from the transmission of infections. In addition to the impermeability requirements of the European standard, the gloves also fulfil the virus impermeability requirements of ASTM F 1671.
- With the newly developed material, the Sempermed[®] Nitrile Skin² feels soft and light.
- The elastic material offers highest wearing comfort and a maximum **sense of touch.**
- The structured fingertips guarantee a sure grip at all times, even in a moist environment.

TECHNICAL DATA

- The glove is powder-free and does not contain any natural latex proteins. This makes it ideal for people with type I allergy who do not want to miss the comfort of latex.
- The Sempermed[®] Nitrile Skin2 is both labelled as medical device (class 1) and single-use protective glove (PPE category III). Thereby it can be used for additional applications such as laboratory/pharma, food/catering, professional cleaning/ hygiene.
- The violet colour makes it easy to distinguish from natural latex gloves and allows for easy recognition.



sempermed[®] ()^{nitrile} skin²

ТҮРЕ	Non-sterile, powder-free examination glove for single use
MATERIAL	Nitrile butadiene rubber (NBR)
COLOUR	Violet
GLOVE SHAPE	Flat with rolled rim
OVERALL LENGHT/SIZE as per EN 455-2	XS, S, M, L, XL median 240 mm
WALL THICKNESS measured in double layer in the palm are	Median 0.14 mm
IMPERABILITY as per EN 455-1	AQL 1.5
FORCE AT BREAK typical production value	Median 6 N
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years
PACKAGING	In dispenser carton XS, S, M, L 100/200 pcs., XL 90/180 pcs. 10 dispenser cartons per shipping carton
MARKING	EN 1041, EN 980, EN 420, ISO 15223
PURPOSE as set out in MDD 93/42 EEC as set out in PPE 89/686 EEC	Medical device class I Singel-use protective glove – PPE category III (protection against chemical substances for limited time) ²⁾

The technical data are average values from the production and may vary in individual cases. ¹⁾ according to EN 455-1 and ASTM F-1671 ²⁾ Please note the Sempermed chemical resistance list, available on www.sempermed.com/sempermed_informs(/chemical_resistance

IMPORTANT NOTE: Actual product information is available under www.sempermed.com. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. Semperit does not accept any liability for incorrect use of the gloves. In case of doubt, obtain expert advice before use. The information and classification correspond to the latest status prior to printing. Technical details are average values from production and may vary in individual cases. Subject to mistakes, printing errors and amendments.



Semperit Technische Produkte Gesellschaft m. b. H. Segment Sempermed Austria · 1030 Vienna · Modecenterstraße 22 Tel.: +43 1 79 777-0 e-mail: sempermed@semperitgroup.com www.sempermed.com